



HydraVer 2 hydrazine reagent, 500 mL

Ürün Numarası:

179049

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For hydrazine determination in boiler waters by the p-Dimethylaminobenzaldehyde method. Hach method 8141. Range: 4-600 µg/L N₂H₄. Also used in visual test kits. Range: 0.02-1.00 mg/L. 500 mL bottle

Teknik özellikler

Cihaz:	All except DR900, PC II
Hacim / Paket boyutu:	500 mL
Metot:	p-Dimethylaminobenzaldehyde
Metot numarası:	8141
Ölçüm aralığı:	4 - 600 µg/L N ₂ H ₄
Parametre:	Hydrazine
Platform:	Çözelti
Saklama koşulları:	10 - 25 °C
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)