



Reagent powder pillows, HR silica, molybdate, 10 mL, pk/1000

Ürün Numarası:

2107328

TRY Fiyat:

Bize Ulaşın

Sevk tarihi gösterilmedi

For determination of high range silica by the silicomolybdate method (Hach method 8185) Range:1 - 100 mg/L. Sample Size: 10 mL sample. Pack of 1000 Powder Pillows.

Teknik özellikler

Hacim / Paket boyutu:	1000 /pk
Metot:	Silicomolybdate
Platform:	Powder Pillows
Test sayısı:	1000

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)