



Reagent powder pillows, silica, citric acid, 10 mL, pk/1000

Ürün Numarası:

2106228

Tehlikeli

TRY Fiyat:

Bize Ulaşın

Sevk tarihi gösterilmedi



For determination of high range silica by the silicomolybdate method (Hach method 8185), and low range silica by the heteropoly blue method (Hach method 8186). Pack of 1000 Powder Pillows.

Teknik özellikler

Hacim / Paket boyutu:	pkg 1000
Parametre:	Silica
Platform:	Powder Pillows

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)