



Reagent solution, amino acid F, 20 ampoules (ultra-low range)

Ürün Numarası:

2386420

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of Ultra Low Range silica by the Heteropoly Blue method. Hach method 8282 Range: 3 - 1000 µg/L SiO₂. Ampuled for increased stability

Teknik özellikler

Metot:	Heteropoly Blue
Platform:	Çözelti
Test sayısı:	20

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)