



## Standard solution, detergent, 60 mg/L as LAS, pk/16 - 10 mL Voluette ampoules

Ürün Numarası:

1427110

TRY Fiyat:  
Mevcut

Bize Ulaşın

For accuracy checks in detergents (surfactants) determinations

### Simple Accuracy Checks

Voluette and PourRite Ampule Standards are designed to work directly with Standard Additions accuracy checks in Hach procedures and photometer software. Use the TenSette Pipet to add precise spikes of standard and check percent recovery, to confirm analytical system performance.

### Minimize waste

Open a fresh standard ampule for each accuracy check. Each ampule contains enough standard for one set of standard additions.

### Reliable results

Sealed ampoules prevent accidental contamination or evaporation of standard solutions.

---

### Teknik özellikler

|                |                             |
|----------------|-----------------------------|
| Konsantrasyon: | 60 mg/L                     |
| Platform:      | Çözelti                     |
| Tanım:         | Detergent Standard Solution |
| Test sayısı:   | 16                          |