



Standard solution, sulfite, 15 mg/L as SO₃ (equivalent), 500 mL

Ürün Numarası:

2408449

TRY Fiyat:
Mevcut

Bize Ulaşın

Sulfite equivalent standard solution. Prepared with stabilized Sodium Thiosulfate. Equivalent to 15 mg/L as SO₃. 500 mL Bottle

Teknik özellikler

Konsantrasyon:

15 mg/L

Platform:

Çözelti

Tanım:

Sulfite Equivalent Solution as SO₃²⁻ (NIST)