



Quality control standard inorganics for drinking water, 500 mL

Ürün Numarası:

2833049

TRY Fiyat:
Mevcut

Bize Ulaşın

Inorganics Quality Control Standard for Drinking Water. Ready to use. No dilution required. Mixed standard includes: Fluoride: 1 mg/L F^- , Nitrate: 2 mg/L NO_3^- -N, Phosphate: 2 mg/L PO_4^{3-} , Sulfate: 50 mg/L SO_4^{2-} . 500 mL Bottle.

Teknik özellikler

Konsantrasyon:

Nitrate 2 mg/L NO_3^- -N

Platform:

Çözelti

Tanım:

Drinking Water Inorganics, Mixed Parameter QC Standard