



Quality control standard hardness, low range, 500 mL

Ürün Numarası: 2833449
TRY Fiyat: Bize Ulaşın
Mevcut

Ready to use - no dilutions required. Mixed standard includes: Total Hardness: 100 mg/L as CaCO₃. Calcium Hardness: 50 mg/L as CaCO₃.

Teknik özellikler

Konsantrasyon: Total Hardness 100 mg/L CaCO₃
Platform: Çözelti
Tanım: Low Range Hardness, Mixed Parameter QC Standard