



Quality control standard metals for drinking water, low range, 500 mL

Ürün Numarası:

2833749

TRY Fiyat:
Mevcut

Bize Ulaşın

Ready to use - no dilutions required. Mixed standard includes: Copper: 1 mg/L Cu, Iron: 0.3 mg/L Fe, Manganese: 0.1 mg/L mg/lMn. 500 mL Bottle.

Teknik özellikler

Konsantrasyon:

Demir: 0,3 mg/L Fe

Platform:

Çözelti

Tanım:

Metals Drinking Water (Low Range) Mixed Parameter QC Standard