



SpecCheck Secondary Gel Standard Kit - Fluoride, 0-2.0 mg/L F

Ürün Numarası:

2712500

TRY Fiyat:
Mevcut

Bize Ulaşın

Solid, ready for use Gel-Standards for instrument verification in the lab and in the field.

Set contains four vials filled with gels that simulate the test colour at various concentrations.

Teknik özellikler

Cihaz:	All except DR/3000, DR 900 and DR/800 Series
Konsantrasyon:	0 - 2,00 mg/L F
Metot:	SPADNS