



Reagent powder pillows, ChromaVer 3, chromium, 5 or 10 mL, pk/1000

Ürün Numarası:

1271028

Tehlikeli

TRY Fiyat:

Bize Ulaşın

Sevk tarihi gösterilmedi



For determination of hexavalent chromium by the Diphenylcarbohydrazide method. USEPA accepted for reporting for wastewater analyses using Hach Method 8023, Diphenylcarbohydrazide Method. Photometric Sample Size: 10 mL. Photometric Range: up to 0.70 mg/L. Kit Sample Size: 5 mL. Kit Range: Up to 1.50 Pack of 1000 powder pillows.

Teknik özellikler

Cihaz:	All except DR/3000, DR/2000, DR 900
EPA uygunluğu:	Evet
Metot:	1,5-Diphenylcarbohydrazide
Metot numarası:	8023
Ölçüm aralığı:	0,010 - 0,700 mg/L Cr ⁶⁺
Parametre:	Chromium
Platform:	Powder Pillows
Test sayısı:	1000

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)