



Reagent powder pillows, CuVer 1, copper, 5 mL, pk/100

Ürün Numarası:

1421199

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of copper by the bicinchoninate method. For use with Hach visual test kits. Range: up to 5.0 mg/L Sample Size: 5 mL

Teknik özellikler

Metot:	Bichinchoninate
Platform:	Powder Pillows

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)