



Reagent powder pillows, hydrosulfite, 25 mL, pk/100

Ürün Numarası:

2118869

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of chelated copper by the bicinchoninate method. Sample Size: 25-mL sample.

Teknik özellikler

Metot:	Bicinchoninate
Platform:	Powder Pillows
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)