



Reagent solution, molybdenum 2, 50 mL MDB

Ürün Numarası:

2352512

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of low range molybdenum by the ternary complex method (Hach method 8169). Range: 0.02-3.00 mg/L Mo 50 mL Marked Dropping Bottle

Teknik özellikler

Hacim / Paket boyutu:	50
Metot:	Ternary Complex
Platform:	Çözelti
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)