



## Reagent powder pillows, NitraVer 5, nitrate, 5 mL, pk/50

Ürün Numarası:

1403546

TRY Fiyat:  
Mevcut

Bize Ulaşın

Tehlikeli



For high range nitrate determination by the NitraVer 5 cadmium reduction method. Range: up to 50 mg/L as nitrate-N. Sample Size: 5 mL. Pack of 50 powder pillows.

### Teknik özellikler

Metot:	Cadmium Reduction
Platform:	Powder Pillows

### Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)