



Reagent powder pillows, NitraVer 6, nitrate, 30 mL, pk/1000

Ürün Numarası: 1411928

Tehlikeli

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.



For determination of low range nitrate by the cadmium reduction method. Range: up to 0.50 mg/L as nitrate-N. Use with NitriVer 3. Sample Size: 30 mL.

Teknik özellikler

Metot: Cadmium Reduction
Platform: Powder Pillows

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)