



Standard solution, sulfuric acid, 19.2 N, 500 mL

Ürün Numarası:

203849

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



Multi-purpose sulfuric acid standard solution. Required for determination of volatile acids by the esterification method. 500 mL bottle

Teknik özellikler

Hacim:	500 mL
Konsantrasyon:	19.2 N
Metot:	Esterification
Platform:	Çözelti

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)