



Reagent set, water-in-oil test 'n tube, 50 tests

Ürün Numarası:

2236945

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For use with WO-1 Water-in-Oil Test Kit, Cat. No. 2237300. Contains reagent vials and sample syringes for 50 tests.

Teknik özellikler

Hacim / Paket boyutu: 25 /pk

Test sayısı: 25

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)