



Reagent set, chloride, LR, digital titrator

Ürün Numarası:

2444100

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat

Tehlikeli



For determination of Chloride by argentometric Digital Titrator titration. Range: 10 - 100 mg/L Cl⁻. Reagent Set contains Chloride 2 Indicator Powder Pillows and Silver Nitrate Digital Titrator Cartridges (2) 0.2256 N. Approximately 100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot:	Silver Nitrate Titration
Metot numarası:	8207
Ölçüm aralığı:	10 - 100 mg/L Cl ⁻
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)