



Sulfide 1 reagent, 500 mL

Ürün Numarası:

181649

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For sulfide determination by the Methylene Blue method (Hach method 8131). Range: up to 800 µg/L. Also used with visual test kits. 500 mL bottle

Teknik özellikler

Metot:	Methylene Blue
Platform:	Çözelti

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)