



Standard solution, potassium iodide-iodate, 0.0125 N, 1 L

Ürün Numarası:

1400153

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For sulfite determination by buret iodometric titration. 0.0125 M (APHA formulation). 1 liter bottle

Teknik özellikler

Metot:	Iodate-Iodide
Platform:	Çözelti
Test sayısı:	40 - 100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)