



Reagent set, sulfite, SU-5 kit

Ürün Numarası:

2440000

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



Reagent set for determination of Sulfite by iodometric drop count titration. Replacement reagent set for SU-5 Test Kit (148002). Range: Up to 200 mg/L SO₃. Set includes Sulfamic Acid Powder Pillows, Sulfite 1 Reagent Powder Pillows and Sulfite 3 Solution for approximately 100 tests.

Teknik özellikler

Önemli özellik:

Sulfite Reagent set, SU-5 Kit : Iodometric Method

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)