



Reagent set, sulfite, buret titration

Ürün Numarası:

2459800

Tehlikeli

TRY Fiyat:

Bize Ulaşın

3 gün içinde teslimat



For determination of Sulfite by iodometric buret titration. Hach Method 8071 Range: 0 to 500+ mg/L SO_3^{2-} . Reagent set includes: Dissolved Oxygen 3 Reagent Powder Pillows, Potassium Iodide-Iodate Standard Solution and Starch Indicator Solution for approximately 100 tests.

Teknik özellikler

Metot numarası:

8071

Ölçüm aralığı:

0 - > 500 mg/L SO_3^{2-}

Önemli özellik:

Sulfite Reagent set, Buret Titration

Test sayısı:

100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)