



## Buffer solution, ULR chlorine, pk/20

Ürün Numarası:

2493120

Tehlikeli

TRY Fiyat:  
Mevcut

Bize Ulaşın



For ultra-low Total Chlorine Determination by the DPD method. Hach Method 8370. Range: 2 - 500 µg/L Cl<sub>2</sub> Ampuled for stability. Package of 20 1.5 mL ampules.

### Teknik özellikler

|                       |         |
|-----------------------|---------|
| Hacim / Paket boyutu: | 20 /pk  |
| Metot:                | DPD     |
| Metot numarası:       | 8370    |
| Platform:             | Çözelti |
| Test sayısı:          | 20      |

### Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)