



Blanking reagent, ULR chlorine, 29 mL

Ürün Numarası:

2493023

Tehlikeli

TRY Fiyat:
Mevcut

Bize Ulaşın



For determination of ultra-low range chlorine by the DPD ULR Method. Range 2 - 500 µg/L, USEPA accepted for reporting wastewater analysis results.

Teknik özellikler

Metot:	DPD
Platform:	Çözelti
Test sayısı:	20

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)