



Reagent, DPD total chlorine, test 'n tube vials, pk/50

Ürün Numarası:

2105645

Tehlikeli

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.



For total chlorine determination by the DPD Method. Test 'N Tube Reagent Vials contain pre-measured DPD reagent. Range: 0.09 - 5.00 mg/L Cl₂. Sample Size: 10 mL

Pack of 50 vials.

Test 'N Tube products provide safe, convenient testing

Capped 16-mm vials provide a self-contained package for mixing and measurement.

Easy ordering

All necessary reagents and vials are contained in the package.

Boric Acid free!

Teknik özellikler

Cihaz:	All except PC II
Metot:	DPD
Metot numarası:	10101
Ölçüm aralığı:	0,09 - 5,00 mg/L Cl ₂
Parametre:	Chlorine, Total
Platform:	Powder Pillows
Test sayısı:	50

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)