



Standard solution, sodium thiosulfate, stabilized, 0.00246 N, 100 mL MDB

Ürün Numarası:

2408532

TRY Fiyat:
Mevcut

Bize Ulaşın

For determination of Chlorine by drop count titration. Used in Hach CN-65 Test Kit (225401). 100 mL Marked Dropping Bottle

Teknik özellikler

Metot:	Iodometric Titration
Platform:	Çözelti

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)