



Reagent set (HR), hypochlorite (bleach), digital titrator

Ürün Numarası:

2687000

TRY Fiyat:

Bize Ulaşın

Sevk tarihi gösterilmedi

Tehlikeli



For determination of Hypochlorite (Bleach) by Digital Titrator titration. Hach Method 10100. Range: 5 - 15% Cl₂ (50 - 150 g/L Cl₂). Stabilized sodium thiosulfate formulation.

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Teknik özellikler

Metot numarası:	10100
Ölçüm aralığı:	5 - 15 % Cl ₂
Platform:	Dijital Titratör
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)