



Reagent set, phosphorus digestion

Ürün Numarası:

2459000

TRY Fiyat:
Mevcut

Bize Ulaşın

Tehlikeli



For acid persulfate digestion of samples for Total Phosphorus (Phosphate) determination by Hach photometric methods. Requires additional reagents for measurement of the reactive phosphorus generated by the digestion. Set contains Potassium Persulfate Powder Pillows, Sodium Hydroxide Solution (5.0 N), and Sulfuric Acid (5.25 N). Approximately 100 tests.

Teknik özellikler

Cihaz:	All
Metot numarası:	8190
Ölçüm aralığı:	0.06-3.50 mg/L
Platform:	Reagent Set
Sindirim gereklidir:	Yes
Test sayısı:	100

Gerekli Aksesuarlar

- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)