



QL filters, pore size 8 - 12 μ m, diameter 110 mm, circular, 100/pk, qualitative

Ürün Numarası:	50656
TRY Fiyat:	Bize Ulaşın
Mevcut	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

Teknik özellikler

Çap:	110 mm
Hacim / Paket boyutu:	100 /pk
Malzeme:	Cellulose