



## Filter, glass fiber, pore size 2.7 µm, diameter 70 mm, 100/pk

|                |             |
|----------------|-------------|
| Ürün Numarası: | 2164300     |
| TRY Fiyat:     | Bize Ulaşın |
| Mevcut         |             |

Made from 100% borosilicate glass fibers. Best used in general laboratory work and as a prefilter for membrane filtration.

### Teknik özellikler

|                       |                  |
|-----------------------|------------------|
| Çap:                  | 70 mm            |
| Gözenek boyutu:       | 2.7 µm           |
| Hacim / Paket boyutu: | 100 /pk          |
| Malzeme:              | Glass Microfiber |