



## Filter, extra thick, 25 mm diameter, OriFlo filter holder, pk/100

Ürün Numarası:

2676700

TRY Fiyat:  
Mevcut

Bize Ulaşın

Designed for rapid, qualitative filtration of particulate matter from aqueous samples using a 25 mm glass fiber filter or other medium with porosity  $\geq 1.0 \mu\text{m}$  (micron). Used in Hach Method 10014 for Ultra Low Range Chlorine in Wastewater (EPA Accepted method).

### Teknik özellikler

|                       |         |
|-----------------------|---------|
| Çap:                  | 25 mm   |
| Hacim / Paket boyutu: | 100 /pk |