



Syringe filter, non-sterile PDVF, 33 mm, 0.45 um, 50/pk

Ürün Numarası:	2513603
TRY Fiyat:	Bize Ulaşın
Mevcut	

Hach syringe filters offer an efficient way to remove particulates from samples prior to analysis. Filters feature 20% more surface area than standard 25 mm filters, which increases flow rates and capacity. PVDF (polyvinylidene fluoride) filters are ideal for solutions containing proteins. Standard Female Luer-Lok inlet with standard male slip Luer outlet. Hydrophilic. Diameter: 33mm. Pore Size: 0.45um. 50/PK.

Teknik özellikler

Birim:	50/pk
Çap:	33 mm
Gözenek boyutu:	0.45 µm
Hacim / Paket boyutu:	50 /pk
Malzeme:	PVDF