



## QL filters, pore size 8 - 12 $\mu$ m, diameter 150 mm, circular, 100/pk, qualitative

Ürün Numarası:	50658
TRY Fiyat:	Bize Ulaşın
Mevcut	

QL Filters are recommended for use in various environmental and wet chemistry reference methods. Extremely pure cellulose filters of nearly 100% alpha content are specially processed to yield negligible extractable levels in the determination and identification of unknown sample materials.

### Teknik özellikler

Çap:	150 mm
Hacim / Paket boyutu:	100 /pk
Malzeme:	Cellulose