



SA filter, cone, diameter 110 mm, soil analysis grade, 500/pk

Ürün Numarası:	2263356
TRY Fiyat:	Bize Ulaşın
Mevcut	

Designed for sandy soils, SA filters are conically folded and have been optimized for use with the common extraction solutions used in soil nutrient analysis (Mehlich 1, Mehlich 3, Morgan, Modified Morgan, Bray, and Olsen). SA Filters provide the best combination of retention rating, filtration speed and purity requirements for soil nutrient testing (P, K, Ca, Mg). Medium filter retention.

Teknik özellikler

Çap:	110 mm
Hacim / Paket boyutu:	50 /pk
Malzeme:	Cellulose