



Flask, Erlenmeyer, glass w/screw cap, 500 mL

Ürün Numarası: 2089649

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Teknik özellikler

Birim:	1 each
Kapasite:	500 mL
Malzeme:	Glass Screw Cap