



## Flask, Erlenmeyer, glass w/screw cap, 1000 mL

Ürün Numarası: 2089653  
TRY Fiyat: Bize Ulaşın  
1 hafta içinde teslimat

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

### Teknik özellikler

Birim: 1 each  
Kapasite: 1000 mL  
Malzeme: Glass Screw Cap