



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 125 mL

Ürün Numarası:

2089673

TRY Fiyat:

Bize Ulaşın

1 hafta içinde teslimat

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Teknik özellikler

Birim:

6/pk

Kapasite:

125 mL

Malzeme:

Glass Screw Cap