



**Be Right™**



## Flask, Erlenmeyer, glass w/screw cap, 6/pk, 250 mL

Ürün Numarası: 2089676

**ÜRETİMDEN KALDIRILAN ÜRÜN**  
Bu ürün artık mevcut değil.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

### Teknik özellikler

Birim:	6/pk
Kapasite:	250 mL
Malzeme:	Glass Screw Cap