



Be Right™



Flask, Erlenmeyer, glass w/screw cap, 6/pk, 1000 mL

Ürün Numarası: 2089683

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.

Graduated glass flask with plastic screw cap. Excellent choice for media preparation, mixing solutions, and culturing. Inert rubber liner cemented into cap and completely autoclavable. Durable graduations indicate approximate volume.

Teknik özellikler

Birim:	6/pk
Kapasite:	1000 mL
Malzeme:	Glass Screw Cap