



Flask, volumetric, polypropylene, 250 mL

Ürün Numarası: 1406046

ÜRETİMDEN KALDIRILAN ÜRÜN
Bu ürün artık mevcut değil.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Teknik özellikler

Birim:	1 each
Gövde malzemesi:	Polypropylene
Kapasite:	250 mL
Sertifikalar:	Not Applicable