



## Beaker, polypropylene, low form, 50 mL

Ürün Numarası: **108041**  
TRY Fiyat: Bize Ulaşın  
Mevcut

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes  $\geq 250$  mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

### Teknik özellikler

Birim:	1 each
Çap:	51.1 mm with spout
Gövde malzemesi:	Polypropylene
Hacim / Paket boyutu:	1 each
Kapasite:	50 mL
Malzeme:	Polypropylene, Low Form