



Beaker, polypropylene, low form, 250 mL

Ürün Numarası: 108046
TRY Fiyat: Bize Ulaşın
Mevcut

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Teknik özellikler

Birim:	1 each
Çap:	83.6 mm with spout
Gövde malzemesi:	Polypropylene
Kapasite:	250 mL
Malzeme:	Polypropylene, Low Form