



## Flask, Erlenmeyer, polymethylpentene capacity 125 mL

Ürün Numarası: 2089843  
TRY Fiyat: Bize Ulaşın  
Mevcut

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

### Teknik özellikler

Birim: 1 each  
Kapasite: 125 mL  
Malzeme: Polymethylpentene  
Tür: Erlenmeyer