



Flask, Erlenmeyer, polymethylpentene capacity 250 mL, 4/pk

Ürün Numarası: 2089876
TRY Fiyat: Bize Ulaşın
1 hafta içinde teslimat

Polymethylpentene flask with polypropylene screwcap features clarity, chemical resistance, strength, and autoclavability (Do not engage threads during cycle). Useful for titrations, preparation and storage of media.

Teknik özellikler

Birim:	4/pk
Kapasite:	250 mL
Malzeme:	Polymethylpentene
Tür:	Erlenmeyer