



## Sample container, 120 mL, pk/250

Ürün Numarası: 2552212  
TRY Fiyat: Bize Ulaşın  
Mevcut

Plastic, leak-proof. Corning Snap Seal offers a non-sterile, leak-proof, attached cap container for field use in sampling water, soils, and sludge. Locking device helps prevent accidental opening of sample; container closes with patented sealing mechanism so sample won't leak.

### Teknik özellikler

Çap: 45 mm (approx)  
Hacim / Paket boyutu: 250 /pk  
Kapasite: 120 mL - 250/pk