



Cylinder, graduated, mixing, glass, 25 mL ± 0.3 mL, 0.5 mL divisions, glass stopper #13

Ürün Numarası:	189640
TRY Fiyat:	Bize Ulaşın
Mevcut	

Low form, single metric scale, glass stopper, graduated. Calibrated "to contain," durable white scale.

Scale is durable white ceramic enamel. Designed to meet ASTM Specification E1272, Style 2 (Class B). Includes glass stopper.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	25 mL
Ölçüm metodu:	To Contain
Sertifikalar:	Class B
Tür:	Tek Metrik, Karıştırma