



Cylinder, graduated, 5 mL $\pm$ 0.1 mL, 0.1 mL divisions (white markings)

Ürün Numarası: 50837
TRY Fiyat: Bize Ulaşın
Mevcut

Calibrated "to deliver" glass cylinder. Single metric scale in durable white markings. Bumper guards included for all sizes $\geq$ 25 mL.

Teknik özellikler

Gövde malzemesi:	Glass
Hacim / Paket boyutu:	1 each
Kapasite:	5 mL
Ölçüm metodu:	To Deliver
Sertifikalar:	Class B
Tür:	Tek Metrik
Yükseklik:	125 mm