



Be Right™



Beaker, griffin, low form, glass, 400 mL

Ürün Numarası:	50048
TRY Fiyat:	Bize Ulaşın
Mevcut	

Griffin style (low form) beakers are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologic

Teknik özellikler

Birim:	1 each
Gövde malzemesi:	Glass
Kapasite:	400 mL
Malzeme:	Glass, Griffin, Low Form