



Be Right™



Bottle, sampling, polycarbonate, 125 mL, 6/pk

Ürün Numarası: 2324333
TRY Fiyat: Bize Ulaşın
Mevcut

Transparent polycarbonate bottle with cap features wide mouth (3.8cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Teknik özellikler

Body Material:	Plastic
Hacim / Paket boyutu:	6 /pk
Kapasite:	125 mL
Malzeme:	PC
Package size UOM:	6pk